

For Research Use Only

LIMPII Recombinant antibody

Catalog Number: 84193-4-RR

Featured Product



Basic Information

Catalog Number:

84193-4-RR

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25908

GenBank Accession Number:

BC021892

GeneID (NCBI):

950

UNIPROT ID:

Q14108

Full Name:

scavenger receptor class B, member 2

Calculated MW:

478 aa, 54 kDa

Observed MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

241468F9

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells, Y79 cells, SH-SY5Y cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84193-4-RR (LIMPII antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



WB result of LIMPII antibody (84193-4-RR; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-LIMPII transfected HepG2 cells.



Biolayer interferometry (BLI) kinetic assays of 84193-4-RR against Human LIMPII were performed. The affinity constant is below 1 pM.