

For Research Use Only

HNRNPAB Recombinant antibody, PBS Only (Detector)

Catalog Number: 84207-3-PBS



Basic Information

Catalog Number:

84207-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35010

GenBank Accession Number:

BC036708

GeneID (NCBI):

3182

UNIPROT ID:

Q99729

Full Name:

heterogeneous nuclear
ribonucleoprotein A/B

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

241462E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01139-2, HNRNPAB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84207-1-PBS. Detection antibody: 84207-3-PBS. Standard: Ag35010. Range: 0.625-80 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84207-3-PBS against Human HNRNPAB was performed. The affinity constant is 16.0 nM.