

For Research Use Only

MLH1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84208-1-PBS



Basic Information

Catalog Number:

84208-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2319

GenBank Accession Number:

BC006850

GeneID (NCBI):

4292

UNIPROT ID:

P40692

Full Name:

mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli)

Calculated MW:

756 aa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

241472C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01140-1, MLH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84208-2-PBS. Detection antibody: 84208-1-PBS. Standard: Ag2319. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01140-3, MLH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84208-3-PBS. Detection antibody: 84208-1-PBS. Standard: Ag2319. Range: 0.313-40 ng/mL



Biolayer Interferometry (BLI) kinetic assay of 84208-1-PBS against Human MLH1 was performed. The affinity constant is 0.642 nM.