

For Research Use Only

S100A4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84221-2-PBS



Basic Information

Catalog Number:

84221-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9019

GenBank Accession Number:

BC016300

GeneID (NCBI):

6275

ENSEMBL Gene ID:

ENSG00000196154

UNIPROT ID:

P26447

Full Name:

S100 calcium binding protein A4

Calculated MW:

101 aa, 12 kDa

Purification Method:

Protein A purification

CloneNo.:

241448B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01155-1, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-3-PBS. Detection antibody: 84221-2-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01155-2, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-2-PBS. Detection antibody: 84221-1-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01155-3, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-2-PBS. Detection antibody: 84221-4-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL