

For Research Use Only

NETO2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84222-1-PBS



Basic Information

Catalog Number: 84222-1-PBS	GenBank Accession Number: BC012381	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81831	CloneNo.: 241451E2
Source: Rabbit	UNIPROT ID: Q8NC67	
Isotype: IgG	Full Name: neuropilin (NRP) and tolloid (TLL)-like	
Immunogen Catalog Number: AG2971	2	
	Calculated MW: 525 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01159-1, NETO2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84222-3-PBS. Detection antibody: 84222-1-PBS. Standard: Ag2971. Range: 0.156-20 ng/mL