

For Research Use Only

NETO2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84222-4-PBS



Basic Information

Catalog Number:

84222-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2971

GenBank Accession Number:

BC012381

GeneID (NCBI):

81831

UNIPROT ID:

Q8NC67

Full Name:

neuropilin (NRP) and tolloid (TLL)-like

2

Calculated MW:

525 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

241451F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01159-2, NETO2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84222-2-PBS. Detection antibody: 84222-4-PBS. Standard: Ag2971. Range: 0.156-20 ng/mL