

For Research Use Only

VDAC2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84225-1-PBS



Basic Information

Catalog Number:	84225-1-PBS	GenBank Accession Number:	BC000165	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	7417	CloneNo.:	241569E1
Source:	Rabbit	UNIPROT ID:	P45880		
Isotype:	IgG	Full Name:	voltage-dependent anion channel 2		
Immunogen Catalog Number:	AG2266	Calculated MW:	294 aa, 32 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01162-1, VDAC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84225-2-PBS. Detection antibody: 84225-1-PBS. Standard: Ag2266. Range: 1.563-200 ng/mL.