

For Research Use Only

NPEPPS Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84230-3-PBS



Basic Information

Catalog Number:

84230-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36035

GenBank Accession Number:

BC065294

GeneID (NCBI):

9520

UNIPROT ID:

P55786

Full Name:

aminopeptidase puromycin sensitive

Purification Method:

Protein A purification

CloneNo.:

241502H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01126-1, NPEPPS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84230-3-PBS. Detection antibody: 84230-1-PBS. Standard: Ag36035. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01126-2, NPEPPS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84230-2-PBS. Detection antibody: 84230-3-PBS. Standard: Ag36035. Range: 0.625-80 ng/mL



Biolayer Interferometry (BLI) kinetic assay of 84230-3-PBS against Human NPEPPS was performed. The affinity constant is 0.364 nM.