

For Research Use Only

TNIK Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84231-2-PBS



Basic Information

Catalog Number:

84231-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34132

GenBank Accession Number:

BC150256

GeneID (NCBI):

23043

UNIPROT ID:

Q9UKE5

Full Name:

TRAF2 and NCK interacting kinase

Purification Method:

Protein A purification

CloneNo.:

241518C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01127-1, TNK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84231-1-PBS. Detection antibody: 84231-2-PBS. Standard: Ag34132. Range: 1.25-80 ng/mL.



Cytometric bead array standard curve of MP01127-2, TNK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84231-2-PBS. Detection antibody: 84231-3-PBS. Standard: Ag34132. Range: 0.781-100 ng/mL.