

For Research Use Only

# Carbonic Anhydrase IX/CA9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84233-6-PBS



## Basic Information

Catalog Number: 84233-6-PBS	GenBank Accession Number: BC014950	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 768	CloneNo.: 241573G12
Source: Rabbit	UNIPROT ID: Q16790	
Isotype: IgG	Full Name: carbonic anhydrase IX	
Immunogen Catalog Number: AG36941	Calculated MW: 459 aa, 50 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01136-3, CA9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84233-1-PBS. Detection antibody: 84233-6-PBS. Standard: Ag36941. Range: 0.313-40 ng/mL