

For Research Use Only

COMMD9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84251-2-PBS



Basic Information

Catalog Number:

84251-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35078

GenBank Accession Number:

BC010892

GeneID (NCBI):

29099

UNIPROT ID:

Q9P000

Full Name:

COMM domain containing 9

Purification Method:

Protein A purification

CloneNo.:

241505B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01110-1, COMMD9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84251-2-PBS. Detection antibody: 84251-3-PBS. Standard: Ag35078. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01110-2, COMMD9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84251-2-PBS. Detection antibody: 84251-1-PBS. Standard: Ag35078. Range: 0.625-40 ng/mL