

For Research Use Only

C12orf11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84252-6-PBS



Basic Information

Catalog Number:

84252-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13727

GenBank Accession Number:

BC003081

GeneID (NCBI):

55726

UNIPROT ID:

Q9NVM9

Full Name:

chromosome 12 open reading frame

11

Calculated MW:

706 aa, 80 kDa

Purification Method:

Protein A purification

CloneNo.:

241586H6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01142-3, ASUN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84252-2-PBS. Detection antibody: 84252-6-PBS. Standard: Ag13727. Range: 0.781-100 ng/mL.