

For Research Use Only

METTL2B Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84257-3-PBS



Basic Information

Catalog Number:

84257-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10410

GenBank Accession Number:

BC107586

GeneID (NCBI):

55798

UNIPROT ID:

Q6P1Q9

Full Name:

methyltransferase like 2B

Calculated MW:

313 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

241490D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01163-2, METTL2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84257-3-PBS. Detection antibody: 84257-2-PBS. Standard: Ag10410. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01163-3, METTL2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84257-1-PBS. Detection antibody: 84257-3-PBS. Standard: Ag10410. Range: 0.781-100 ng/mL