

For Research Use Only

# GBP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84264-4-PBS



## Basic Information

Catalog Number:

84264-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2425

GenBank Accession Number:

BC022272

GeneID (NCBI):

2634

UNIPROT ID:

P32456

Full Name:

guanylate binding protein 2,  
interferon-inducible

Calculated MW:

591 aa, 67 kDa

Purification Method:

Protein A purification

CloneNo.:

241580E12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01177-2, GBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84264-4-PBS. Detection antibody: 84264-3-PBS. Standard: Ag2425. Range: 0.156-20 ng/mL