

For Research Use Only

Rat IL-1 beta Recombinant antibody, PBS Only (Capture)

Catalog Number: 84274-2-PBS



Basic Information

Catalog Number:

84274-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031512

GeneID (NCBI):

24494

UNIPROT ID:

Q63264

Full Name:

interleukin 1 beta

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

241001G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01210-1, RAT IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84274-2-PBS. Detection antibody: 84274-3-PBS. Standard: Eg31371. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01210-2, RAT IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84274-2-PBS. Detection antibody: 84274-4-PBS. Standard: Eg31371. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01210-3, RAT IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84274-2-PBS. Detection antibody: 84274-1-PBS. Standard: Eg31371. Range: 0.313-40 ng/mL