

For Research Use Only

# MLPH Recombinant antibody, PBS Only (Capture)

Catalog Number: 84282-3-PBS



## Basic Information

Catalog Number:

84282-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32931

GenBank Accession Number:

BC001653

GeneID (NCBI):

79083

UNIPROT ID:

Q9BV36

Full Name:

melanophilin

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

241621H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01196-2, MLPH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84282-3-PBS. Detection antibody: 84282-5-PBS. Standard: Ag32931. Range: 0.313-40 ng/mL