

For Research Use Only

NOG Recombinant antibody, PBS Only (Capture)

Catalog Number: 84283-6-PBS



Basic Information

Catalog Number:

84283-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32755

GenBank Accession Number:

BC034027

GeneID (NCBI):

9241

UNIPROT ID:

Q13253

Full Name:

noggin

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

241511E5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01174-4, Human Noggin Recombinant Matched Antibody Pair - PBS only. 84283-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag32755. 84283-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL



The mean Noggin concentration was determined to be 79.96 ng/mL in mouse brain tissue extract based on a 2.0 mg/mL extract load.