

For Research Use Only

# PYROXD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84289-5-PBS



## Basic Information

Catalog Number:

84289-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34768

GenBank Accession Number:

BC021662

GeneID (NCBI):

79912

UNIPROT ID:

Q8WU10

Full Name:

pyridine nucleotide-disulphide  
oxidoreductase domain 1

Purification Method:

Protein A purification

CloneNo.:

241639G7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01201-3, PYROXD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84289-4-PBS. Detection antibody: 84289-5-PBS. Standard: Ag34768. Range: 0.313-40 ng/mL.