

For Research Use Only

ADAM17 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84292-1-PBS



Basic Information

Catalog Number:

84292-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32418

GenBank Accession Number:

BC136783

GeneID (NCBI):

6868

UNIPROT ID:

P78536

Full Name:

ADAM metallopeptidase domain 17

Calculated MW:

824 aa, 93 kDa

Purification Method:

Protein A purification

CloneNo.:

241305D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01169-1, ADAM17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84292-2-PBS. Detection antibody: 84292-1-PBS. Standard: Ag32418. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01169-2, ADAM17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84292-3-PBS. Detection antibody: 84292-1-PBS. Standard: Ag32418. Range: 0.625-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84292-1-PBS against Human ADAM17 was performed. The affinity constant is 12.8 nM.