

For Research Use Only

ARFGEF2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84296-2-PBS



Basic Information

Catalog Number:

84296-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35396

GenBank Accession Number:

BC050449

GeneID (NCBI):

10564

UNIPROT ID:

Q9Y6D5

Full Name:

ADP-ribosylation factor guanine
nucleotide-exchange factor 2
(brefeldin A-inhibited)

Purification Method:

Protein A purification

CloneNo.:

241538F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01199-1, ARFGEF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84296-2-PBS. Detection antibody: 84296-3-PBS. Standard: Ag35396. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01199-2, ARFGEF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84296-2-PBS. Detection antibody: 84296-1-PBS. Standard: Ag35396. Range: 0.313-40 ng/mL