

For Research Use Only

# NMNAT2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84307-1-PBS



## Basic Information

Catalog Number:

84307-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26909

GenBank Accession Number:

BC020998

GeneID (NCBI):

23057

UNIPROT ID:

Q9BZQ4

Full Name:

nicotinamide nucleotide  
adenylyltransferase 2

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

241484B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01207-2, NMNAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84307-1-PBS. Detection antibody: 84307-2-PBS. Standard: Ag26909. Range: 0.156-20 ng/mL.