

For Research Use Only

Mouse 4-1BB/TNFRSF9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84315-1-PBS



Basic Information

Catalog Number:

84315-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001077509.1

GeneID (NCBI):

21942

UNIPROT ID:

P20334

Full Name:

tumor necrosis factor receptor
superfamily, member 9

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

241647D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01216-2, MOUSE 4-1BB/TNFRSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84315-1-PBS. Detection antibody: 84315-3-PBS. Standard: Eg1627. Range: 0.078-10 ng/mL