

For Research Use Only

# HNF4A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84316-4-PBS



## Basic Information

Catalog Number: 84316-4-PBS	GenBank Accession Number: BC137539	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3172	CloneNo.: 241593C8
Source: Rabbit	UNIPROT ID: P41235	
Isotype: IgG	Full Name: hepatocyte nuclear factor 4, alpha	
Immunogen Catalog Number: AG24647	Calculated MW: 53 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01183-2, HNF4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84316-4-PBS. Detection antibody: 84316-5-PBS. Standard: Ag24647. Range: 0.781-100 ng/mL.