

For Research Use Only

ATPBD3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84318-2-PBS



Basic Information

Catalog Number:	84318-2-PBS	GenBank Accession Number:	NM_145232	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	90353	CloneNo.:	241630E9
Source:	Rabbit	UNIPROT ID:	Q7Z7A3		
Isotype:	IgG	Full Name:	ATP binding domain 3		
Immunogen Catalog Number:	AG35037	Calculated MW:	36 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01189-1, ATPBD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84318-2-PBS. Detection antibody: 84318-3-PBS. Standard: Ag35037 Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01189-2, ATPBD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84318-2-PBS. Detection antibody: 84318-1-PBS. Standard: Ag35037. Range: 0.625-80 ng/mL.