

For Research Use Only

A1CF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84322-2-PBS



Basic Information

Catalog Number:	84322-2-PBS	GenBank Accession Number:	BC054873	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	29974	CloneNo.:	241644C9
Source:	Rabbit	UNIPROT ID:	Q9NQ94		
Isotype:	IgG	Full Name:	APOBEC1 complementation factor		
Immunogen Catalog Number:	AG22929	Calculated MW:	594 aa, 65 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01198-1, A1CF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84322-2-PBS. Detection antibody: 84322-3-PBS. Standard: Ag22929. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01198-3, A1CF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84322-2-PBS. Detection antibody: 84322-1-PBS. Standard: Ag22929. Range: 0.313-40 ng/mL.