

For Research Use Only

UTY Recombinant antibody, PBS Only (Detector)

Catalog Number: 84332-2-PBS



Basic Information

Catalog Number:

84332-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35806

GenBank Accession Number:

NM_007125.4

GeneID (NCBI):

7404

UNIPROT ID:

O14607

Full Name:

ubiquitously transcribed
tetra-tryptophan repeat gene, Y-
linked

Calculated MW:

150kDa, 1347aa

Purification Method:

Protein A purification

CloneNo.:

241626A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01230-3, UTY Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84332-3-PBS. Detection antibody: 84332-2-PBS. Standard: Ag35806. Range: 0.313-40 ng/mL

Cytometric bead array standard curve of MP01230-1, UTY Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84332-1-PBS. Detection antibody: 84332-2-PBS. Standard: Ag35806. Range: 0.313-40 ng/mL