

For Research Use Only

# UTY Recombinant antibody, PBS Only (Capture)

Catalog Number: 84332-4-PBS



## Basic Information

Catalog Number:	84332-4-PBS	GenBank Accession Number:	NM_007125.4	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	7404	CloneNo.:	241626D6
Source:	Rabbit	UNIPROT ID:	O14607		
Isotype:	IgG	Full Name:	ubiquitously transcribed tetratricopeptide repeat gene, Y-linked		
Immunogen Catalog Number:	AG35806	Calculated MW:	150kDa, 1347aa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01230-2, UTY Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84332-4-PBS. Detection antibody: 84332-3-PBS. Standard: Ag35806. Range: 0.313-40 ng/mL