

For Research Use Only

C5AR1/CD88 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84336-1-PBS



Basic Information

Catalog Number: 84336-1-PBS	GenBank Accession Number: NM_001736.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 728	CloneNo.: 241712H5
Source: Rabbit	UNIPROT ID: P21730	
Isotype: IgG	Full Name: complement component 5a receptor 1	
	Calculated MW: 39kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01255-1, C5AR1/CD88 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84336-4-PBS. Detection antibody: 84336-1-PBS. Standard: Eg2185. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01255-3, C5AR1/CD88 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84336-2-PBS. Detection antibody: 84336-1-PBS. Standard: Eg2185. Range: 0.313-40 ng/mL.