

For Research Use Only

INTS1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84339-1-PBS



Basic Information

Catalog Number:

84339-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35428

GenBank Accession Number:

BC069262

GeneID (NCBI):

26173

UNIPROT ID:

Q8N201

Full Name:

integrator complex subunit 1

Purification Method:

Protein A purification

CloneNo.:

241627F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01233-1, INTS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84339-1-PBS. Detection antibody: 84339-3-PBS. Standard: Ag35428. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01233-2, INTS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84339-1-PBS. Detection antibody: 84339-2-PBS. Standard: Ag35428. Range: 0.625-80 ng/mL