

For Research Use Only

ORAI3 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84340-3-PBS



Basic Information

Catalog Number:

84340-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36135

GenBank Accession Number:

BC006126

GeneID (NCBI):

93129

UNIPROT ID:

Q9BRQ5

Full Name:

ORAI calcium release-activated
calcium modulator 3

Calculated MW:

295 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

241568A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01228-1, ORAI3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84340-3-PBS. Detection antibody: 84340-2-PBS. Standard: Ag36135. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01228-2, ORAI3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84340-1-PBS. Detection antibody: 84340-3-PBS. Standard: Ag36135. Range: 0.625-80 ng/mL