

For Research Use Only

GAPVD1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84343-2-PBS



Basic Information

Catalog Number:	84343-2-PBS	GenBank Accession Number:	BC021119	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	26130	CloneNo.:	241609C3
Source:	Rabbit	UNIPROT ID:	Q14C86		
Isotype:	IgG	Full Name:	GTPase activating protein and VPS9 domains 1		
Immunogen Catalog Number:	AG10492	Calculated MW:	1478 aa, 165 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01218-1, GAPVD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84343-2-PBS. Detection antibody: 84343-1-PBS. Standard: Ag10492. Range: 0.781-100 ng/mL.