

For Research Use Only

Placental lactogen Recombinant antibody, PBS Only (Detector)

Catalog Number: 84361-1-PBS



Basic Information

Catalog Number:

84361-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0201

GenBank Accession Number:

BC002717

GeneID (NCBI):

1442

UNIPROT ID:

PODML3

Full Name:

chorionic somatomammotropin hormone 1 (placental lactogen)

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

241625C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01252-1, CSH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84361-2-PBS. Detection antibody: 84361-1-PBS. Standard: Ag0201. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01252-3, CSH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84361-3-PBS. Detection antibody: 84361-1-PBS. Standard: Ag0201. Range: 0.781-100 ng/mL.