

For Research Use Only

Mouse SIGNR1/CD209b Recombinant antibody, PBS Only (Capture)

Catalog Number: 84365-2-PBS



Basic Information

Catalog Number:

84365-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_026972.5

GeneID (NCBI):

69165

UNIPROT ID:

Q8CJ91-1

Full Name:

CD209b antigen

Calculated MW:

37kDa

Purification Method:

Protein A purification

CloneNo.:

241726C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01231-1, MOUSE SIGNR1/CD209b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84365-2-PBS. Detection antibody: 84365-3-PBS. Standard: Eg1497. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01231-2, MOUSE SIGNR1/ CD209b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84365-2-PBS. Detection antibody: 84365-1-PBS. Standard: Eg1497. Range: 0.625-80 ng/mL