

For Research Use Only

POLR3C Recombinant antibody, PBS Only (Detector)

Catalog Number: 84367-3-PBS



Basic Information

Catalog Number:

84367-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27311

GenBank Accession Number:

BC002586

GeneID (NCBI):

10623

UNIPROT ID:

Q9BUI4

Full Name:

polymerase (RNA) III (DNA directed)

polypeptide C (62kD)

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

241402C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01234-1, POLR3C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84367-4-PBS. Detection antibody: 84367-3-PBS. Standard: Ag27311. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01234-3, POLR3C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84367-1-PBS. Detection antibody: 84367-3-PBS. Standard: Ag27311. Range: 0.781-100 ng/mL