

For Research Use Only

POLR3C Recombinant antibody, PBS Only (Capture)

Catalog Number: 84367-5-PBS



Basic Information

Catalog Number:

84367-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27311

GenBank Accession Number:

BC002586

GeneID (NCBI):

10623

UNIPROT ID:

Q9BUI4

Full Name:

polymerase (RNA) III (DNA directed)

polypeptide C (62kD)

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

241402F12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01234-4, Human POLR3C Recombinant Matched Antibody Pair - PBS only. 84367-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27311. 84367-4-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean POLR3C concentration was determined to be 6.1 ng/mL in HEK-293T cell extract based on a 1.4 mg/mL extract load and 8.9 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load.