

For Research Use Only

# POLR3C Recombinant antibody, PBS Only

Catalog Number: 84367-7-PBS



## Basic Information

Catalog Number:

84367-7-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27311

GenBank Accession Number:

BC002586

GeneID (NCBI):

10623

UNIPROT ID:

Q9BUI4

Full Name:

polymerase (RNA) III (DNA directed)

polypeptide C (62kD)

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

241402D3

## Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using POLR3C antibody (84367-7-RR, Clone: 241402D3 ) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84367-7-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84367-7-RR against Human POLR3C were performed. The affinity constant is below 1 pM.