

For Research Use Only

# Mouse CD28 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84368-5-PBS



## Basic Information

Catalog Number:

84368-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_007642.4

GeneID (NCBI):

12487

UNIPROT ID:

P31041

Full Name:

CD28 antigen

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

241707H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01248-1, MOUSE CD28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84368-5-PBS. Detection antibody: 84368-3-PBS. Standard: Eg1408. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01248-3, MOUSE CD28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84368-5-PBS. Detection antibody: 84368-4-PBS. Standard: Eg1408. Range: 0.313-40 ng/mL.