

For Research Use Only

# Rat CD44 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84369-1-PBS



## Basic Information

Catalog Number:

84369-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

XM\_006234631.3

GeneID (NCBI):

25406

UNIPROT ID:

O70509

Full Name:

Cd44 molecule

Calculated MW:

40kDa

Purification Method:

Protein A purification

CloneNo.:

241648B6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP01236-1, RAT CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84369-2-PBS. Detection antibody: 84369-1-PBS. Standard: Eg1458. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01236-2, RAT CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84369-3-PBS. Detection antibody: 84369-1-PBS. Standard: Eg1458. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01236-3, RAT CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84369-4-PBS. Detection antibody: 84369-1-PBS. Standard: Eg1458. Range: 0.156-20 ng/mL