

For Research Use Only

THAP11 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84373-1-PBS



Basic Information

Catalog Number:	84373-1-PBS	GenBank Accession Number:	BC012182	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	57215	CloneNo.:	241590A5
Source:	Rabbit	UNIPROT ID:	Q96EK4		
Isotype:	IgG	Full Name:	THAP domain containing 11		
Immunogen Catalog Number:	AG19342	Calculated MW:	314 aa, 34 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP01244-3, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-1-PBS. Detection antibody: 84373-2-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP01244-3, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-1-PBS. Detection antibody: 84373-2-PBS. Standard: Ag19342. Range: 0.781-100 ng/mL