

For Research Use Only

THAP11 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84373-4-PBS



Basic Information

Catalog Number:

84373-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19342

GenBank Accession Number:

BC012182

GeneID (NCBI):

57215

UNIPROT ID:

Q96EK4

Full Name:

THAP domain containing 11

Calculated MW:

314 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

241590C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01244-2, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-4-PBS. Detection antibody: 84373-3-PBS. Standard: Ag19342. Range: 0.781-100 ng/mL.