

For Research Use Only

HAPLN1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84374-3-PBS



Basic Information

Catalog Number:

84374-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25221

GenBank Accession Number:

BC057808

GeneID (NCBI):

1404

UNIPROT ID:

P10915

Full Name:

hyaluronan and proteoglycan link protein 1

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

241588F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01245-2, HAPLN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84374-3-PBS. Detection antibody: 84374-1-PBS. Standard: Ag25221. Range: 0.313-40 ng/mL.