

For Research Use Only

# COX18 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84377-2-PBS



## Basic Information

Catalog Number:

84377-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12709

GenBank Accession Number:

BC101684

GeneID (NCBI):

285521

UNIPROT ID:

Q8N8Q8

Full Name:

COX18 cytochrome c oxidase assembly homolog (S. cerevisiae)

Calculated MW:

333 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

241359E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01250-1, COX18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84377-2-PBS. Detection antibody: 84377-1-PBS. Standard: Ag12709. Range: 0.781-100 ng/mL.