

For Research Use Only

# MYH13 Recombinant antibody

Catalog Number: 84379-2-RR



## Basic Information

Catalog Number:

84379-2-RR

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_003802

GeneID (NCBI):

8735

UNIPROT ID:

Q9UIK3

Full Name:

myosin, heavy chain 13, skeletal muscle

Calculated MW:

224 kDa

Observed MW:

230-240 kDa

Purification Method:

Protein A purification

CloneNo.:

241179A3

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse eye tissue, rat eye tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84379-2-RR (MYH13 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.