

For Research Use Only

INTS8 Recombinant monoclonal antibody, PBS Only

Catalog Number: 84382-2-PBS



Basic Information

Catalog Number:

84382-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34627

GenBank Accession Number:

BC136754

GeneID (NCBI):

55656

UNIPROT ID:

Q75QN2

Full Name:

integrator complex subunit 8

Calculated MW:

111 kDa

Observed MW:

111 kDa

Purification Method:

Protein A purification

CloneNo.:

253187D5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Integrator complex subunit 8 (INTS8) is one of the major components of RNA polymerase II and has been demonstrated to be involved in the cleavage of small nuclear RNAs and transcriptional processes (PMID: 34880328). INTS8 was robustly increased in neurodevelopmental diseases and numerous tumors. High INTS8 expression has a tight correlation with poor prognosis in multi-cancers (PMID: 30881114). Western blot analysis detected INTS8 at an apparent molecular mass of 111 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84382-2-RR (INTS8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84382-2-PBS in a different storage buffer formulation.