

For Research Use Only

PDGFR alpha/CD140a Recombinant antibody, PBS Only (Detector)

Catalog Number: 84383-3-PBS



Basic Information

Catalog Number:

84383-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015186

GeneID (NCBI):

5156

UNIPROT ID:

P16234

Full Name:

platelet-derived growth factor
receptor, alpha polypeptide

Calculated MW:

123 kDa

Purification Method:

Protein A purification

CloneNo.:

241658H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01256-2, PDGFR alpha/CD140a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84383-2-PBS. Detection antibody: 84383-3-PBS. Standard: RP00003. Range: 0.781-100 ng/mL