

For Research Use Only

# PDGFRA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84383-5-PBS



## Basic Information

Catalog Number:

84383-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015186

GeneID (NCBI):

5156

UNIPROT ID:

P16234

Full Name:

platelet-derived growth factor  
receptor, alpha polypeptide

Calculated MW:

123 kDa

Purification Method:

Protein A purification

CloneNo.:

241659E9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01256-1, PDGFR alpha/CD140a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84383-5-PBS. Detection antibody: 98141-1-PBS. Standard: RP00003. Range: 0.781-100 ng/mL