

For Research Use Only

# TF Recombinant antibody, PBS Only (Detector)

Catalog Number: 84386-10-PBS



## Basic Information

Catalog Number:

84386-10-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1317

GenBank Accession Number:

BC020671

GeneID (NCBI):

7018

UNIPROT ID:

P02787

Full Name:

transferrin

Purification Method:

Protein A purification

CloneNo.:

250421B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Transferrin/TF concentration was determined to be 1.9 mg/mL in human serum.



Cytometric bead array standard curve of MP01258-4, Transferrin/TF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84386-11-PBS. Detection antibody: 84386-10-PBS. Standard: HZ-1317. Range: 0.156-40 ng/mL.