

For Research Use Only

Transferrin/TF Recombinant antibody, PBS Only (Detector)

Catalog Number: 84386-2-PBS



Basic Information

Catalog Number:

84386-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC059367

GeneID (NCBI):

7018

UNIPROT ID:

P02787

Full Name:

transferrin

Calculated MW:

698 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

241769A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01258-1, Transferrin/TF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84386-4-PBS. Detection antibody: 84386-2-PBS. Standard: Eg0114. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01258-3, Transferrin/TF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84386-3-PBS. Detection antibody: 84386-2-PBS. Standard: Eg0114. Range: 0.781-100 ng/mL.