

For Research Use Only

IL18R1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84390-7-PBS



Basic Information

Catalog Number:

84390-7-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093975

GeneID (NCBI):

8809

UNIPROT ID:

Q13478

Full Name:

interleukin 18 receptor 1

Calculated MW:

541 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

243056G10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01261-4, Human IL-18R1 Recombinant Matched Antibody Pair - PBS only. 84390-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard RP00006. 84390-6-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean IL-18R1 concentration was determined to be 8,042.0 pg/mL in NK-92 cell extract based on a 2.0 mg/mL extract load and 1,797.2 pg/mL in HEK-293 cell extract based on a 2.2 mg/mL extract load.