

For Research Use Only

IL-12RB1/CD212 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84391-1-PBS



Basic Information

Catalog Number:

84391-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC029121

GeneID (NCBI):

3594

UNIPROT ID:

P42701

Full Name:

interleukin 12 receptor, beta 1

Calculated MW:

662 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

241660B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01262-3, IL12RB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84391-3-PBS. Detection antibody: 84391-1-PBS. Standard: RP00005. Range: 0.781-50 ng/mL.